

AmnioWorks™

CellRight Technologies®



AmnioWorks™

AmnioWorks™ is derived from human amniotic membrane recovered following a live birth from volunteer donors. Amnion tissue is known to be an immune-privileged tissue containing growth factors and cytokines that impact tissue repair and regeneration. AmnioWorks™ is available in multiple sizes designed to closely optimize membrane defects.

AmnioWorks™ amnion patch grafts are sterilized using low-dose, low-temperature gamma irradiation to achieve a sterility assurance level (SAL) of 10^{-6} . The grafts should be stored in ambient temperatures and have a shelf life of 5-years. AmnioWorks™ amnion patch may be rehydrated according to the clinician's preference, and may be used as membrane for applications in accordance with the physician's well-informed medical judgment.

To place an order:

customercare@cellrighttechnologies.com

T: 210.659.9353

F: 210.659.9556

© 2018 CellRight Technologies, LLC.

MKT-0030-A 05/2018

