

Conventional Grafts

CellRight Technologies®

As trusted stewards of anatomical donations, CellRight Technologies® makes every effort to maximize this precious gift.



Tricortical Block



Unicortical Block

CellRight Technologies® proprietary BioRinse™ technology extends to the processing of conventional grafts including Unicortical and Tricortical Blocks. These structural grafts can be used in spinal fusion, foot and ankle, and general orthopedic reconstruction surgeries, according to a physician's well-informed medical judgement.


The grafts are freeze-dried and sterilized using a validated low-dose, low-temperature gamma irradiation to achieve a sterility assurance level (SAL) of 10^{-6} . The grafts should be stored in ambient temperatures (59-86°F or 15-30°C) and have a shelf-life of 5 years from the date of packaging.



CellRight Technologies® employs strict quality assurance and quality control procedures to ensure patient safety. CellRight's medical director, a licensed physician, performs an extensive medical review of the donor's medical/social history to determine eligibility. CellRight Technologies® only accepts donors whose screening, serologic, and microbiologic tests meet or exceed the current standards established by the Food and Drug Administration (FDA) and the American Association of Tissue Banks (AATB).


Conventional grafts may be used in surgical procedures including spine, neuro, orthopedics, trauma, foot and ankle, reconstruction, and dental.

| PRODUCT CODE | PRODUCT NAME | PRODUCT DIMENSION |
|--------------|---------------------------|--|
| 421018 | Tricortical Block 16-18mm | (L)≥25mm, (W)=16-18mm, (T1)>5mm,(T2)>5mm |
| 421022 | Tricortical Block 20-22mm | (L)≥25mm, (W)=20-22mm, (T1)>5mm,(T2)>5mm |
| 421026 | Tricortical Block 24-26mm | (L)≥25mm, (W)=24-26mm, (T1)>5mm,(T2)>5mm |



Tricortical Block

| PRODUCT CODE | PRODUCT NAME | PRODUCT DIMENSION |
|--------------|------------------------|-------------------------------------|
| 421005 | Unicortical Block 5mm | (L)=11-13mm, (W) =13-15mm, (T)=5mm |
| 421006 | Unicortical Block 6mm | (L)=11-13mm, (W) =13-15mm, (T)=6mm |
| 421007 | Unicortical Block 7mm | (L)=11-13mm, (W) =13-15mm, (T)=7mm |
| 421008 | Unicortical Block 8mm | (L)=11-13mm, (W) =13-15mm, (T)=8mm |
| 421009 | Unicortical Block 9mm | (L)=11-13mm, (W) =13-15mm, (T)=9mm |
| 421010 | Unicortical Block 10mm | (L)=11-13mm, (W) =13-15mm, (T)=10mm |



Unicortical Block

To place an order and/or to inquire about additional conventional graft sizes please contact customer care@cellrighttechnologies.com or call: 210.659.9353



CELLRIGHT[®]
TECHNOLOGIES

1808 Universal City Blvd.
Universal City, Texas 78148
Tel: 210.659.9353
Fax: 210.659.9556

www.cellrighttechnologies.com